

长兴广达新型材料有限公司 **Changxing Guangda New Materials CO.,LTD**

GDGEL-282 RHEOLOGICAL ADDITIVE

Description

GDGEL-282 is an organo clay (modified bentonite), designed especially for use in solvent based system. It can be used widely in solvents from low polarity to medium-high polarity.

Properties

Composition Organic derivative of a montmorillonite clay

Fine powder, light yellow **Appearance**

Moisture content (105°C,2hr) 3.5% Max Loss on ignition (1000°C) 35% Max

Bulk density 0.35-0.45g/cm3

Fineness (<74 µ m. %) 95% Min

- **Advantages** High gelling efficiency, high viscosity
 - Wide application range, it can be applied in all low, middle and high polarity range
 - Produce reproducible thixotropic consistency over a wide temperature range
 - Imparts particle suspension, preventing hard settling of pigment and fillers
 - Exerts strong film reinforcing action in organic binder systems

Application

- Marine paint, Industrial paint, heavy-duty, anticorrosion paint such as bituminous coatings, polyester paints, Alkyd paints.
- Sealant
- Inks
- Grease
- Oil Drilling Mud

Using

GDGEL-282 disperses under high shear and requires a polar activator for best efficiency. Heat is not required to obtain good dispersion. Addition levels are typically in 0.2%~2.0%. Polar activator (95% Ethanol or 95% Methanol) is needed. Normally addition level of polar activator is recommended to be 45%~60% of GDGEL-282 in weight. If in xylene system, addition level of polar activator is about 30%-40% of GDGEL-282 in weight.

Add: Wushangang Town, Changxing County, Huzhou City, Zhejiang, China

Website: www.gdgel.com E-mail: gdgel@gdgel.com

TEL: 86-572-6513637 FAX:86-572-6513637 Pre-gel is made as below procedure.

- 1. Charge solvent to mill
- 2. Slowly add GDGEL-282 (6-10% by weight of total pre-gel) to mill under agitation.
- 3. Mix at high speed for 5 minutes.
- 4. Add polar activator.
- 5. Mix at high speed for 5 minutes.

Package &Storage

Kraft paper bag with PE liner or multi-wall paper sacks or customized. 25kg/bag or 50LB/bag or customized. GDGEL-282 has a shelf life of 24 months, and should be stored in ventilated cool and dry place.

Notice

The information on use is based on data which are believed reliable, but any recommendation or suggestion made is without guarantee or warranty, since the conditions of use are outside our control. All products are sold on the conditions that purchasers shall make their own tests to determine the suitability of such products for their purpose and that all risks are assumed by user. We disclaim any responsibility for damages resulting from careless or improper handling or use. Nothing herein is to be taken as permission, inducement or recommendation to practice any patented invention without a license.

Add: Wushangang Town, Changxing County, Huzhou City, Zhejiang, China

Website: www.gdgel.com E-mail: gdgel@gdgel.com TEL: 86-572-6513637 FAX:86-572-6513637